

Autoclave



Stainless steel steam sterilizer with built-in pressure and temperature safety for hospitals, clinics, labs and industry.

Safety systems and key features

- Automatic shutdown system for overtemperature and overpressure
- Safety relief valve for pressure protection
- High-grade stainless steel sterilization chamber
- Uniform heat distribution inside the chamber
- Heavy-duty structure designed for continuous operation

Applications

- Hospitals and surgical centers
- Medical and research laboratories
- Pharmaceutical manufacturing facilities
- Scientific research institutions
- Industrial sterilization processes requiring high precision

Models and specifications

Model	Capacity (L)	Control system	Sterilization temperature (°C)	Operating pressure	Cycle time	Power (kW)	Power supply
CLX20	20	Manual	126	0.14–0.16 MPa	15–30 min	2	220V / 50 Hz
CLX35	35	Digital	126	0.14–0.16 MPa	20–40 min	3	220V / 50 Hz
CLX75	75	Digital	137	2.5 bar	25–50 min	4.5	220V / 50 Hz



Model	Capacity (L)	Control system	Sterilization temperature (°C)	Operating pressure	Cycle time	Power (kW)	Power supply
CLX100	100	Digital	137	2.5 bar	30–60 min	6	220V / 50 Hz

Warranty

One-year comprehensive warranty on all devices · Two-year warranty against manufacturing defects · Maintenance contract after the warranty period