

Laboratory cabinet system



Modular laboratory cabinet system with floor and wall cabinets, drawer modules, a continuous worktop and an integrated sink.

Key features

- Modular configuration allowing units to be arranged according to space and operational requirements
- Durable black worktop offering a spacious and easy-to-clean working surface
- Integrated sink and mixer tap for convenient washing and laboratory preparation activities
- Multiple storage options, including double-door cabinets and four-drawer units
- Wall-mounted cabinets provide additional storage while keeping the worktop clear
- Ergonomic layout for easy access to tools, materials, and equipment
- Elevated plinth base improves stability and facilitates floor cleaning
- Strong metal handles and fittings designed for frequent daily use
- Bright orange finish creates a modern, distinctive, and highly visible appearance
- Clean, streamlined design suitable for professional laboratory and healthcare environments

Worktop materials

The bench includes a chemical-resistant worktop, usually made from epoxy resin, phenolic resin, compact laminate or black granite, depending on the required chemical and heat resistance.

Warranty

Five-year warranty against manufacturing defects from the acceptance date