

DOUBLE WATER DISTILLATION UNIT

Double-stage distillation for ultra-pure water with maximum efficiency and safety

OVERVIEW

The LabX Double Water Distillation Unit is expertly designed to produce ultra-pure distilled water through a two-stage distillation system, perfectly suited for laboratories, hospitals, and industrial settings.

The unit is made entirely from high-grade stainless steel, with an electrostatically coated outer body for additional protection against humidity and corrosion.

It also features an internal water storage tank to collect the distilled water before being dispensed, ensuring a smooth and consistent flow.



INCLUDED FEATURES:

- 1 Dual Independent Boilers** | Operate in parallel to ensure greater efficiency and reliability.
- 2 Internal Collection Tank** | Stores distilled water before dispensing via external outlet.
- 3 Electrostatically Coated Exterior** | Protects against rust and moisture, easy to clean.
- 4 Full Safety Systems** | Mechanical water sensors, overheat protection, dry operation shut-off, and pressure relief valve.
- 5 Efficient Cooling System** | Built-in heat exchanger reduces cooling water consumption.

LabXTech
QUALITY YOU CAN TRUST