



FLASH POINT TESTER

OVERVIEW

The LabX Flash Point Tester is designed to determine the flash point of liquids safely and accurately. It is widely used in laboratories and industries to assess the ignition temperature of flammable substances, ensuring compliance with safety standards.



INCLUDED FEATURES:

6 DIGITAL DISPLAY
Provides precise temperature readings

5 EASY OPERATION
Simple controls for user-friendly operation

4 SAFETY
Built-in safety mechanisms to prevent accidents during testing

3 IGNITION SOURCE
Gas ignition system for reliable flame detection

2 TEMPERATURE RANGE
Measures up to 350°C

1 OPEN CUP METHOD
Uses the open cup technique for flash point measurement