

# SINGLE WATER DISTILLATION UNIT

Precision, Reliability, and Efficiency in Every Drop

## OVERVIEW

The LabX Water Distillation Unit is the ideal solution for obtaining pure, impurity-free water suitable for scientific laboratories, medical facilities, and industrial applications.

Constructed entirely from 316L stainless steel, and equipped with advanced protection systems, the unit ensures long-lasting performance, high purity output, and safe operation.



## INCLUDED FEATURES:

- 1 High-Quality Construction** Built entirely of stainless steel to resist corrosion and ensure hygiene.
- 2 Advanced Protection System** Automatic safety shutoff for the heating element in case of low water level.
- 3 Automatic Water Level Control** Maintains continuous operation without manual intervention.
- 4 Reliable Performance** Consistent and efficient water distillation using advanced technology

**LabXTech**  
QUALITY YOU CAN TRUST