



WATER BATH

Precision, Reliability, and Efficiency in Every Drop

OVERVIEW

The LabX Tech Water Bath is a precise laboratory instrument designed to heat and maintain samples at constant temperatures using water as the heat transfer medium.

It is widely used in laboratories for chemical, biological, pharmaceutical, and educational applications.

With digital temperature control, high-quality construction, and multiple safety features, this unit ensures excellent performance and durability.



INCLUDED FEATURES:

- 1 Microprocessor PID temperature controller
- 2 Built-in temperature calibration and digital display
- 3 Timer function for auto shut-off
- 4 Float sensor for water level monitoring
- 5 Stainless steel bath within a sturdy body
- 6 Electrostatic powder-coated body (anti-corrosion)
- 7 Over-temperature protection and auto-tuning
- 8 Novel structure prevents vapor leakage
- 9 Compact, industrial-grade build
- 10 Easy to use, clean, and maintain
- 11 Long service life
- 12 Uniform heating distribution