



# BOD INCUBATOR

## OVERVIEW

The BOD Incubator from LabX Tech offers reliable and accurate temperature-controlled conditions for performing BOD (Biochemical Oxygen Demand) testing. It's specifically engineered to meet the rigorous standards of environmental laboratories, water treatment facilities, academic institutions, and research organizations.



## INCLUDED FEATURES:

7

### ECO-FRIENDLY INSULATION:



Polyurethane foam for thermal stability and energy efficiency.

6

### ACCURATE TEMPERATURE CONTROL:



PID digital controller for temperature range 5°C to 60°C.

5

### SMART ALARM SYSTEM:



Visual and audio alerts for temperature deviations or faults.

4

### LED DISPLAY:



Shows real-time internal temperature.

2

### UNIFORM AIRFLOW SYSTEM:



Internal fan ensures consistent heat distribution.

1

### HYGIENIC & DURABLE INTERIOR:



304 stainless steel chamber with 3 to 5 adjustable shelves

3

### MOBILITY:



Strong castor wheels for easy transport.